

HRK-MATE SHORT INTRODUCTION



Think Customer & Market Group

Contents

1. **About HRK-Mate**
2. **Connectivity in HRK-Mate**
3. **Preview**

1.


About HRK-Mate

Operator



HRK Series
HNT-1/1P

Transmit the data to Physician's PC



Physician




HRK-Mate

- Doctor can review the measurement data obtained from the HRK Series and HNT-1/1P through the HRK-Mate software.

2.

Connectivity in HRK-Mate



3.

3.1. HRK Series

-. REF & KER Information

Patient Information									
Pat. NO	P00033			Exam Date	2013-02-13				
Name	Kim	Jungshan			Pat. ID	kim598			Edit
Gender	☒ Male ☐ Female			Birth Date	1983-11-19			Delete	
Note	It's me								


REF & KER Information									
[REF] VD : 12.0, PD : 24					[KER] Index : 1.3375				
<R>	SPH	CYL	AX		<L>	SPH	CYL	AX	
	-2.98	-1.39	89			-2.87	-1.41	90	
	-2.98	-1.39	89			-2.87	-1.41	90	
	-2.98	-1.40	89			-2.87	-1.41	90	
AVG	-2.98	-1.39	89		AVG	-2.87	-1.41	90	

<R>	R1	R2	AX		<L>	R1	R2	AX	
	7.94	7.78	93			7.94	7.77	91	
	7.94	7.78	93			7.94	7.77	91	
	7.94	7.78	93			7.94	7.77	91	
	mm	D	AX			mm	D	AX	
R1	7.94	42.50	93		R1	7.94	42.50	91	
R2	7.78	43.40	3		R2	7.77	43.45	1	
AVG	7.86	42.95			AVG	7.86	42.97		
CYL		-0.90	93		CYL		-0.95	91	

3. Preview

3.1. HRK Series

-. Wave Front (Zernike map and coefficients)



3. Preview

3.1. HRK Series


-. PSF(Point Spread Function) & Retinal Image



3. Preview

3.1. HRK Series


- Illumination Mode and Color View Image



3. Preview

3.1. HRK Series




- . Meibo TBUT Image



3. Preview

3.2. HNT-1/1P




- . Pachymetry Image



3. Preview

3.2. HNT-1/1P

-. ACA Image





THANK YOU